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MEDOP EARPLUGS SIGILO PLUS / SIGILO PLUS DETECTABLE EN

MEDOP earplugs SIGILO PLUS comply with Standard EN 352-2:2002. They are reusable. These earplugs has been manufactured according the 352-2:2002 standard, preventing the ear from becoming completely occluded, and the worker from becoming isolated from the environment, whilst at the same time the worker cannot hear any buzzing as a result of the occlusion. The earplugs are manufactured using TPR with a nominal diameter between 8 and 12 mm, which stop humidity from settling in the earplugs, thus eliminating any possible appearance of micro-organisms and irritations of the ear canal. These materials do not cause any diseases, allergic reactions or any other harmful effects on the skin. On the contrary, they provide great softness and comfort. The cord joining the two earplugs enables them to be carried around the neck when the wearer is outside areas of high noise.

1. USAGE

In order to obtain maximum protection, the earplug must be placed correctly. It is recommended that the earplugs be cleaned after each use with soap and water. These earplugs should not be used in areas where the cord could be hooked during the use.

Check the earplugs and replace them if necessary:

- Grasp the cylindrical part of the earplug.
- With the other hand pressed over your head, pull the top of the ear up and out in order to straighten the ear duct.
- Gently put the earplug in place with a light twisting movement.
- To remove the earplug, hold the base and pull out twisting gently at the same time.

2. CARE INSTRUCTIONS

Wash the ear protectors regularly in lukewarm water using neutral soap. Rinse with clean lukewarm water and leave to dry in the air.

Store the earplugs, preferably inside the package, in a ventilated, cool place, avoiding humidity and dust. The earplugs must be inspected regularly, and replaced if damaged or worn. At regular periods, depending on the amount of use, disinfected the unit using a product offering protection against a wide range of bacteria.

NOTE: This equipment is for personal use, and must not be used by any worker other than the owner.

3. ACOUSTIC REDUCTION VALUES

The earplugs should be worn continuously in noise zones. The noise damping levels obtained in official tests were as follows:

SIGILO PLUS

3.1. Noise reduction by Frequency bands in Hz (test complying with ISO/DIS 4869-2)								
Frequency in Hz	63	125	250	500	1.000	2.000	4.000	8.000
Average Reduction	31,6	32,7	30,9	30,6	37,4	29,9	40,3	44,6
Standard Deviation	5,6	5,7	5,9	5,2	8,1	5,1	7,7	5,0
Achieved Reduction	26,0	27,0	25,0	25,4	29,3	24,8	32,6	39,6

3.2. Global reduction in dB

High (H)=27 dB - Medium (M)=27dB - Low (L)=26dB

3.3. Acoustic Reduction Values: SNR=29dB

SIGILO PLUS DETECTABLE

3.1. Noise reduction by Frequency bands in Hz (test complying with ISO/DIS 4869-2)								
Frequency in Hz	63	125	250	500	1.000	2.000	4.000	8.000
Average Reduction	31,4	31,3	29,2	31,4	36,6	30,3	40,8	44,7
Standard Deviation	5,5	4,8	4,8	4,0	6,2	5,7	7,3	5,1
Achieved Reduction	25,9	26,5	24,4	27,4	30,4	24,6	33,5	39,6

3.5. Global reduction in dB

High (H)=27 dB - Medium (M)=28dB - Low (L)=27dB

3.6. Acoustic Reduction Values: SNR=30dB

IMPORTANT: The noise reductions can only be achieved provided that the ear protectors:

- Are placed correctly, adjusted and cared for according to instructions in this manual.
- Are used continuously in noisy areas, and are checked regularly replacing them if damaged.

NOTE: The best ear protector will only be effective if used during the whole duration of the exposure to noise. Removing the protector even for a short period of time, seriously reduces protection. The following recommendations are not followed: The protection provided by these earplugs will be significantly reduced, thus reducing their effectiveness. These earplugs can be adversely affected by certain chemicals. If you wish further information please contact the manufacturer. These earplugs should not be used when there is a risk of hooking during its using.

4. RECOMMENDATIONS

SUITABLE PROTECTION: The best protection would result in an actual exposure of 75 to 80 dBA.

OVERPROTECTION: It is not advisable to try and obtain greater noise reduction than required, as overprotection involves greater discomfort for the user. The protecting effect is better if a merely sufficient noise reduction level is chosen.

WARNING: These earplugs may be adversely affected by certain chemicals. If you wish further information please carry out the activity correctly or for the worker's safety may be obscured.

WARNING: These instructions describe correct use of the equipment, and must be followed. The equipment may be adversely affected by certain chemicals. If you wish further information please contact the manufacturer. These earplugs should not be used when there is a risk of hooking during its using.

NOTE: This equipment is for personal use, and must not be used by any worker other than the owner.

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